

AUTO Refractometer

Matronix Q30+ Auto Refractometer

The most advanced technique of image processing.
Stylish Design, Multi Function Keys, Movement Lock,
The precision optical technique, consistent measurement
and the accurate readings.

Fully Automatic



Technical Features:

Refraction measurement	Sphere(sph) : -25D~+22D、 0.12D/0.25D steps Cylinder(cyl) : -8D~+8D、 0.12D/0.25D steps Axis : 0~180°、 1° steps
Pupil distance	45~85mm、 0.5mm steps
Vertical distance	0m、 12mm、 13.75mm
Minimum pupil diameter	Φ2.0mm
Charter	Automatic fogging method
Extra Features	IOL Mode
Display	5.7inch true color LCD
Print	Impact less thermal printer
Power supply	100~120V/200~240V、 50/60HZ
Dimension	275(W)*475(D)*435-465(H)mm
Weight	18kg



Sales and Service in India by:

BioVision,

**#12, 1st L Main Road, 2nd Stage, 1st Block,
Nagarabhavi, Bangalore-560072**

T: 080-23281009; 080-65371500

M: 9448468510

E-Mail: bv_helpline@hotmail.com

Web: www.biovision.co.in