



Mfg.:- Ophthalmic Chair Unit, Operation Table, Motorized Table, Eye Testing Drum, Surgeon Stool, Surgical Blade & Instrument

OPHTHALMIC CHAIR UNIT MAX-2007

FEATURES

UNIT - Smooth-rotating, Well- balanced Sliding Table Top with ample space for two instruments ~ Automatic, under-the-table connection for slit Lamp ~ Touch-activated Safety Switch to stop inadvertent chair movements ~ Integrated, centralized Control Console with touchstone Keys ~ Strong Pole to support Lamp, Phoropter, Projector & Near Vision Holder ~ Convenient Instrument Holder Tray ~ Dust-free, enclosed Trial Lens Rack with ~ straight, unobtrusive, in & out sliding movement.

CONTROL PANEL

Elegant, single-piece Console panel for centralize unit control. ~ Multifaceted control for chair movements ~ Built-in connectivity for Indirect Ophthalmoscope, Near Vision Holder & Projector.

CHAIR.

Fully flexible, rotating recliner chair ~ *silent, CE-marked, German DC motor-driven, up-down movements and inclination*, and Zero position Memory ~ Adjustable footrest, which can be fully raised to sleeping position ~ Foot Switch for fluid up / down movement control ~ Comfortable polyurethane-filled Seat ~ Full movement control from back rest panel.


An Electro Mechanical Refraction Unit

The Compact Elegant Comprehensive Integration Unit



SPECIFICATIONS

Power Supply AC	230 V AC
Power Consumption	100 W
Load Lifting - Up/Down	135 Kg
Vertical Traveling-of chair	200 mm.
Seat Height minimum	580 mm.
Seat Height maximum	780 mm.
Net Weight-with Pack	240 Kg. (Approx)
Minimum Area Requir	2325 mm. X 1420mm.



Ordering Information - Available Dr/Right side & Left side & Colour



Colour Shown in this Picture Give you an approximate Idea.

Sales & Service in India by:



No. 12, 1st 'L' Main, 2nd Stage 1st Block, Nagarabhavi, Bangalore 560 072.
 Phone : 080 23281009, E-mail: bv_helpline@hotmail.com